



# Needle Holders

## 07

## 07 Needle Holders

Adson	94
Brown	94
Collier	93
Crile	98
Crile-Murray	94
Crile-Wood	94
Derf	93
Finochietto	96
Gillies	97
Hegar-Baumgartner	93
Halsey	93
Heaney	95
Jamison	95
Johnson	96
Lichtenberg	98
M.G. Hospital	94
Masson	95
Mathieu	98
Mayo-Hegar	95
Metzenbaum	94
Needle Holder	93-98
Neivert	93
Olsen-Hegar	97
Stratte	97
Toennis	98
Wangensteen	96
Webster	93

Needle Holders  
Porta-agujas  
Porte-aiguilles



**SS-07-101**  
12 cm

one fenestrated jaw  
un lado de la boca  
fenestrada  
uncote de mors fenetre



Derf



**SS-07-102**  
13 cm

Serration  
boca con estrias  
longitudinales.  
Mors avec rainures  
Longitudinales



Collier



**SS-07-103**  
12,5 cm  
grooved jaw  
SSa acanalada  
mors strie



Webster




**SS-07-104**  
13 cm  
smooth jaw  
boca lisa  
mors lisse



**SS-07-105**  
13,5 cm  
smooth jaw  
boca lisa  
mors lisse



Halsey



**SS-07-106**  
12,5 cm  
smooth jaw, special groove  
boca lisa, especialmente acanalada  
mors lisse, avce stries speciales



Neivert



**SS-07-107**  
14,5 cm



Hegar-Baumgartner



**SAPTESTAR**  
 Manufacturer & Exporter Of Surgical, Dental & Veterinary Instruments

Needle Holders  
 Porta-agujas  
 Porte-aiguilles




**SS-07-108**  
 13 cm

Heavy pattern  
 modelo fuerte  
 modele fort

Diagonally serrated jaw  
 boca estriada transversalmente  
 mors strie transversalment

**Brown**



**SS-07-109**  
 18 cm

**Crile-Murry**



**SS-07-110**  
 15 cm  
 slender pattern  
 modelo delgado  
 modele gracile



**SS-07-111**  
 15 cm

**Crile-Wood**



**SS-07-112**  
 18 cm  
 one fenestrated jaw  
 un lado de la boca fenestrada  
 un cote de mors fenetre

**Adson**



**SS-07-113**  
 18,5 cm  
 one grooved and fenestrated jaw  
 un lado de la boca canalado y fenestrado  
 un cote de mors strie et fenetre


**Metzenbaum**



**SS-07-114**  
 17.5 cm  
 heavy pattern  
 modelo fuerte  
 modele fort

**M.G Hospital**

Needle Holders  
 Porta-agujas  
 Porte-aiguilles





Needle Holders  
Porta-agujas  
Porte-aiguilles



Finochietto



Wangensteen




Johnson

Needle Holders  
 Porta-agujas  
 Porte-aiguilles



**SS-07-131**  
 23 cm

**Stratte**



**SS-07-132**  
 14 cm / 5 1/2"

**SS-07-133**  
 17 cm / 6 1/4"

**Olsen-Hegar**



**SS-07-134**  
 16 cm

**Gillies**

Needle Holders  
 Porta-agujas  
 Porte-aiguilles



Crile



1/4

**SS-07-135**  
 15 cm



Mathieu



**SS-07-136**  
 14 cm



**SS-07-137**  
 17 cm



**SS-07-138**  
 20 cm



Mathieu



**SS-07-139**  
 14 cm



1/4

**SS-07-140**  
 16.5 cm



Toennis



1/4

**SS-07-141**  
 18.5 cm



Lichtenberg



1/4

**SS-07-142**  
 20 cm